

# Human CPB2 (His tag) recombinant protein



Catalog Number: 519543

## General Information

### Gene Name Synonym

Carboxypeptidase B2, Carboxypeptidase U, CPU, Plasma Carboxypeptidase B, pCPB, Thrombin-Activable Fibrinolysis Inhibitor, TAFI, CPB2

### Protein Construction

A DNA sequence encoding the Human CPB2 Phe23-Val423 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 1mM ZnCl<sub>2</sub>, 10% Glycerol, pH7.5.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human CPB2 (His tag) recombinant protein