

# Human SOAT2/ACACT2/ACAT2 (Trx & His tag) recombinant protein



Catalog Number: 519591

## General Information

### Gene Name Synonym

Acetyl-CoA acetyltransferase cytosolic, Acetyl-CoA transferase-like protein, Cytosolic acetoacetyl-CoA thiolase, ACTL, ACAT2

### Protein Construction

A DNA sequence encoding the Human SOAT2/ACACT2/ACAT2 Met1-Glu397 is expressed with a Trx, 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 20mM PB□ 150mM NaCl,pH7.4.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human SOAT2/ACACT2/ACAT2 (Trx & His tag)  
recombinant protein