

Mouse Angptl4/Farp/Fiaf/Ng27 (Fc tag) recombinant protein



Catalog Number: 519633

General Information

Gene Name Synonym

Angiopoietin-related protein
4,425O18-1,Angiopoietin-like protein 4,Fasting-
induced adipose factor,Hepatic
fibrinogen/angiopoietin-related
protein,HFARP,Secreted protein
Bk89,Angptl4,Farp, Fiaf, Ng27

Protein Construction

A DNA sequence encoding the Mouse
Angptl4/Farp/Fiaf/Ng27 Lys167-Ser410 is
expressed with a Fc tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL
method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of
PBS,pH7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Mouse Angptl4/Farp/Fiaf/Ng27 (Fc tag)
recombinant protein