

# Human MMP1/CLG (His tag) recombinant protein



Catalog Number: 519652

## General Information

### Gene Name Synonym

Interstitial Collagenase, Fibroblast Collagenase, Matrix Metalloproteinase-1, MMP-1, MMP1, CLG

### Protein Construction

A DNA sequence encoding the Human MMP1/CLG Phe20-Asn469 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 20mM MES, 150mM NaCl, 0.05% Brij35, pH 5.5.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human MMP1/CLG (His tag) recombinant protein