

# Human CPB1/CPB/PCPB (His tag) recombinant protein



Catalog Number: 519672

## General Information

### Gene Name Synonym

Carboxypeptidase B, Pancreas-Specific Protein, PASP, CPB1, CPB, PCPB

### Protein Construction

A DNA sequence encoding the Human CPB1/CPB/PCPB His16-Tyr417 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu$ g protein as determined by the LAL method.

## Stability

Store for up to 6 months at -20°C

## Formulation


Supplied as a 0.2  $\mu$ m filtered solution of 20mM MES, 150mM NaCl, pH 5.5.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

## SDS-PAGE



Human CPB1/CPB/PCPB (His tag) recombinant protein