

Human CFB/BF/BFD (His tag) recombinant protein



Catalog Number: 519673

General Information

Gene Name Synonym

Complement Factor B, C3/C5 Convertase, Glycine-Rich Beta Glycoprotein, GBG, PBF2, Properdin Factor B, CFB, BF, BFD

Protein Construction

A DNA sequence encoding the Human CFB/BF/BFD Thr26-Leu764 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human CFB/BF/BFD (His tag) recombinant protein