

# Human FKBP3/FKBP25 (His tag) recombinant protein



Catalog Number: 519688

## General Information

### Gene Name Synonym

Peptidyl-prolyl cis-trans isomerase FKBP3, PPIase FKBP3, 25 kDa FK506-binding protein, 25 kDa FKBP, FKBP-25, FK506-binding protein 3, FKBP-3, Immunophilin FKBP25, Rapamycin-selective 25 kDa immunophilin, Rotamase, FKBP25

### Protein Construction

A DNA sequence encoding the Human FKBP3/FKBP25 Met1-Asp224 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu$ g protein as determined by the LAL method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2  $\mu$ m filtered solution of 20mM TrisHCl, 1mM DTT, 10% glycerol.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human FKBP3/FKBP25 (His tag) recombinant protein