

Human FKBP3/FKBP25 (His tag) recombinant protein



Catalog Number: 519688

General Information

Gene Name Synonym

Peptidyl-prolyl cis-trans isomerase FKBP3, PPIase FKBP3, 25 kDa FK506-binding protein, 25 kDa FKBP, FKBP-25, FK506-binding protein 3, FKBP-3, Immunophilin FKBP25, Rapamycin-selective 25 kDa immunophilin, Rotamase, FKBP25

Protein Construction

A DNA sequence encoding the Human FKBP3/FKBP25 Met1-Asp224 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM TrisHCl, 1mM DTT, 10% glycerol.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human FKBP3/FKBP25 (His tag) recombinant protein