

# Human TINAGL1/GIS5/LCN7/OLRG2 (His tag) recombinant protein



Catalog Number: 519692

## General Information

### Gene Name Synonym

Tubulointerstitial nephritis antigen-like (TINAGL1), also known as Glucocorticoid-inducible protein 5, Oxidized LDL-responsive gene 2 protein, TIN Ag-related protein, Tubulointerstitial nephritis antigen-related protein, GIS5, LCN7, OLRG2 and TINAGL. belongs to the peptidase C1 family.

### Protein Construction

A DNA sequence encoding the Human TINAGL1/GIS5/LCN7/OLRG2 Ala22-His467 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu$ g protein as determined by the LAL method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 20mM PB, 150mM NaCl, 1mM EDTA, 2m MDTT, pH7.4.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Human TINAGL1/GIS5/LCN7/OLRG2 (His tag)  
recombinant protein