

# Human NECAP2 (His tag) recombinant protein



Catalog Number: 519707

## General Information

### Gene Name Synonym

Adaptin ear-binding coat-associated protein 2,  
NECAP endocytosis associated protein 2, NECAP2

### Protein Construction

A DNA sequence encoding the Human NECAP2  
Met1-Phe263 is expressed with a 6His tag at the  
N-terminus, 6His tag at the C-terminus.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu$ g protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2  $\mu$ m filtered solution of 20mM PB $\square$   
150mM NaCl,pH7.4.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human NECAP2 (His tag) recombinant protein