

Human FH recombinant protein



Catalog Number: 519709

General Information

Gene Name Synonym

Fumarate Hydratase Mitochondrial, Fumarase, FH

Protein Construction

A DNA sequence encoding the Human FH Ala45-Lys510 is expressed.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of 20mM Tris, pH 8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human FH recombinant protein