

# Human SULT1C2/SULT1C1 (His tag) recombinant protein



Catalog Number: 519716

## General Information

### Gene Name Synonym

Sulfotransferase 1C2, ST1C2, Sulfotransferase 1C1, SULT1C#1, humSULTC2, SULT1C2, SULT1C1,

### Protein Construction

A DNA sequence encoding the Human SULT1C2/SULT1C1 Met1-Leu296 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 20mM Tris, 0.1M NaCl, pH 8.5.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human SULT1C2/SULT1C1 (His tag) recombinant protein