

# Human PTPRN/ICA3/ICA512 (His tag) recombinant protein



Catalog Number: 519722

## General Information

### Gene Name Synonym

Receptor-type tyrosine-protein phosphatase-like N, R-PTP-N, Islet cell antigen 512, ICA 512, Islet cell autoantigen 3, PTP IA-2, PTPRN, ICA3, ICA512

### Protein Construction

A DNA sequence encoding the Human PTPRN/ICA3/ICA512 Gln607-Asn686&Trp795-Leu888 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu$ g protein as determined by the LAL method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2  $\mu$ m filtered solution of PBS, pH 7.4.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human PTPRN/ICA3/ICA512 (His tag)  
recombinant protein