

# Human SH2D1A/DSHP/SAP (His tag) recombinant protein



Catalog Number: 519742

## General Information

### Gene Name Synonym

SH2 Domain-Containing Protein 1A, Duncan Disease SH2-Protein, Signaling Lymphocytic Activation Molecule-Associated Protein, SLAM-Associated Protein, T-Cell Signal Transduction Molecule SAP, SH2D1A, DSHP, SAP

### Protein Construction

A DNA sequence encoding the Human SH2D1A/DSHP/SAP Met1-Pro128 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, 10% Glycerol, pH 7.5.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human SH2D1A/DSHP/SAP (His tag) recombinant protein