

# Human PTPRN/ICA3/ICA512 (His tag) recombinant protein



Catalog Number: 519743

## General Information

### Gene Name Synonym

Receptor-type tyrosine-protein phosphatase-like N, R-PTP-N, Islet cell antigen 512, ICA 512, Islet cell autoantigen 3, PTP IA-2, PTPRN, ICA3, ICA512

### Protein Construction

A DNA sequence encoding the Human PTPRN/ICA3/ICA512 Arg576-Gln950 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL

method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 20mM Tris, 150mM NaCl, pH8.0.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Human PTPRN/ICA3/ICA512 (His tag)  
recombinant protein