

# Human DTD1/C20orf88/DUEB/HARS2 (His tag) recombinant protein



Catalog Number: 519754

## General Information

### Gene Name Synonym

D-tyrosyl-tRNA(Tyr) deacylase 1, DNA-unwinding element-binding protein B, DUE-B, Histidyl-tRNA synthase-related, DTD1, C20orf88, DUEB, HARS2

### Protein Construction

A DNA sequence encoding the Human DTD1/C20orf88/DUEB/HARS2 Met1-Pro209 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human DTD1/C20orf88/DUEB/HARS2 (His tag)  
recombinant protein