

Human KLK6/PRSS18/PRSS9 (His tag) recombinant protein



Catalog Number: 519800

General Information

Gene Name Synonym

Kallikrein-6, Neurosin, Protease M, SP59, Serine Protease 18, Serine Protease 9, Zyme, KLK6, PRSS18, PRSS9

Protein Construction

A DNA sequence encoding the Human KLK6/PRSS18/PRSS9 Glu17-Lys244 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 5mM HCl, 150mM NaCl.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human KLK6/PRSS18/PRSS9 (His tag)
recombinant protein