

# Human LAMTOR3/MAP2K1IP1/MAPKSP1/PRO2783 (His tag) recombinant protein



Catalog Number: 519815

## General Information

### Gene Name Synonym

Ragulator Complex Protein LAMTOR3, Late Endosomal/Lysosomal Adaptor and MAPK and MTOR Activator 3, MEK-Binding Partner 1, Mp1, Mitogen-Activated Protein Kinase Kinase 1-Interacting Protein 1, Mitogen-Activated Protein Kinase Scaffold Protein 1, LAMTOR3, MAP2K1IP1, MAPKSP1

### Protein Construction

A DNA sequence encoding the Human LAMTOR3/MAP2K1IP1/MAPKSP1/PRO2783 Met1-Val124 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, 5mM DTT, 30% Glycerol, pH 8.0.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human  
LAMTOR3/MAP2K1IP1/MAPKSP1/PRO2783 (His  
tag) recombinant protein