

Mouse Col3a1 (His tag) recombinant protein



Catalog Number: 519852

General Information

Gene Name Synonym

Collagen alpha-1(III) chain, Col3a1

Protein Construction

A DNA sequence encoding the Mouse Col3a1 Gln155-Gly1219 is expressed with a 6His tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM HAc-NaAc, 150mM NaCl, pH4.5.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Mouse Col3a1 (His tag) recombinant protein