

Human HSPB7/CVHSP (His tag) recombinant protein



Catalog Number: 519862

General Information

Gene Name Synonym

Pregnancy specific beta-1-glycoprotein 2, PSBG2, Pregnancy-specific beta-1 glycoprotein E, pregnancy-specific beta-1-glycoprotein 2, pregnancy-specific beta-1-glycoprotein 7, PSG2

Protein Construction

A DNA sequence encoding the Human HSPB7/CVHSP Met1-Ile170 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM Tris, 0.2M NaCl, 2mM DTT, 50% glycerol, pH 8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human HSPB7/CVHSP (His tag) recombinant protein