

Human SULT1B1/ST1B2/SULT1B2 (His tag) recombinant protein



Catalog Number: 519872

General Information

Gene Name Synonym

Sulfotransferase Family Cytosolic 1B Member 1, ST1B1, Sulfotransferase 1B1, Sulfotransferase 1B2, ST1B2, Thyroid Hormone Sulfotransferase, SULT1B1, ST1B2, SULT1B2

Protein Construction

A DNA sequence encoding the Human SULT1B1/ST1B2/SULT1B2 Met1-Ile296 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM Tris, 0.1M NaCl, pH 8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human SULT1B1/ST1B2/SULT1B2 (His tag)
recombinant protein