

Human PPID/CYP40/CYPD (His tag) recombinant protein



Catalog Number: 519894

General Information

Gene Name Synonym

Peptidyl-Prolyl Cis-Trans Isomerase D, PPIase D, 40 kDa Peptidyl-Prolyl Cis-Trans Isomerase, Cyclophilin-40, CYP-40, Cyclophilin-Related Protein, Rotamase D, PPID, CYP40, CYPD

Protein Construction

A DNA sequence encoding the Human PPID/CYP40/CYPD Met1-Ala370 is expressed with a 6His tag at the N-terminus, 6His tag at the C-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM PB, 150mM NaCl, pH7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human PPID/CYP40/CYPD (His tag) recombinant protein