

Human RETN/FIZZ3/HXCP1/RSTN (His tag) recombinant protein



Catalog Number: 519915

General Information

Gene Name Synonym

Resistin, Adipose tissue-specific secretory factor, Cysteine-rich secreted protein FIZZ3, C/EBP-epsilon-regulated myeloid-specific secreted cysteine-rich protein, Cysteine-rich secreted protein A12-alpha-like 2, FIZZ3, HXCP1, RSTN, RETN

Protein Construction

A DNA sequence encoding the Human RETN/FIZZ3/HXCP1/RSTN Lys19-Pro108 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted

samples are stable at < -20°C for 3 months.

Formulation

Lyophilized from a 0.2 µm filtered solution of 20mM Acetic acid, pH 3.0.

Usage Guide


Storage

Lyophilized protein should be stored at <-20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days.

Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE



Human RETN/FIZZ3/HXCP1/RSTN (His tag)
recombinant protein