

Human CANX (His tag) recombinant protein



Catalog Number: 519918

General Information

Gene Name Synonym

Calnexin, IP90, Major Histocompatibility Complex Class I Antigen-Binding Protein p88, p90, CANX

Protein Construction

A DNA sequence encoding the Human CANX His21-Pro481 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 2mM CaCl, 10% Glycerol, pH 7.5.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human CANX (His tag) recombinant protein