

# Human NOL3/ARC/NOP (His tag) recombinant protein



Catalog Number: 519923

## General Information

### Gene Name Synonym

Nucleolar Protein 3, Apoptosis Repressor With CARD, Muscle-Enriched Cytoplasmic Protein, Myp, Nucleolar Protein of 30 kDa, Nop30

### Protein Construction

A DNA sequence encoding the Human NOL3/ARC/NOP Met1-Ser208 is expressed.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

## Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

## Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 25mM TrisHCl, 1mM DTT, 1mM EDTA, 2mM  $\beta$ -ME, 20% Glycerol, pH 7.5.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

## SDS-PAGE



Human NOL3/ARC/NOP (His tag) recombinant protein