

# Human NRG1/GGF/HGL/HRGA/NDF/SMDF recombinant protein



Catalog Number: 519957

## General Information

### Gene Name Synonym

Pro-neuregulin-1, Neuregulin-1 beta 1, NRG1-beta 1, HRG1-beta 1, EGF, NRG1, GGF, HGL, HRGA, NDF, SMDF,

### Protein Construction

A DNA sequence encoding the Human NRG1/GGF/HGL/HRGA/NDF/SMDF Thr176-Lys246 is expressed.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Store for up to 12 months at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$  as lyophilized powder. Aliquots of reconstituted samples are stable at  $< -20^{\circ}\text{C}$  for 3 months.

## Formulation

Lyophilized from a  $0.2\ \mu\text{m}$  filtered solution of 20mM PB, 150mM NaCl, pH7.4.

## Usage Guide


### Storage

Lyophilized protein should be stored at  $< -20^{\circ}\text{C}$ , though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at  $4-7^{\circ}\text{C}$  for 2-7 days.

### Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than  $100\ \mu\text{g}/\text{ml}$ . Dissolve the lyophilized protein in ddH<sub>2</sub>O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

## SDS-PAGE



Human NRG1/GGF/HGL/HRGA/NDF/SMDF recombinant protein