

Rat Il2rb (XM_006241970) cDNA/ORF clone



Catalog Number: 708388-2

General Information

Gene Name:

interleukin 2 receptor subunit beta

Official Symbol: Il2rb

Organism: Rattus norvegicus

RefSeq: XM_006241970

Description

Sequence Description:

Identical with the Gene Bank Ref. ID sequence.

Vector: pcDNA3.1-3xFlag

Restriction Sites: HindIII + KpnI

Shipping carrier:

Each tube contains approximately 5 µg - 10 µg of lyophilized plasmid.

Storage:

The lyophilized plasmid can be stored at ambient temperature for three months.

Quality control:

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

Sequencing primer list:

T7:TAATACGACTCACTATAGGG

BGH:TAGAAGGCACAGTCGAGG

Plasmid Resuspension protocol

1. Centrifuge at 5,000×g for 5 min.
2. Carefully open the tube and add 20 µl of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than 5000×g.
5. Store the plasmid at -20 °C.

The plasmid is ready for:

Restriction enzyme digestion; PCR amplification; E. coli transformation; DNA sequencing


E.coli strains for transformation (recommended but not limited):

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5α and TOP10F'.

Vector Information

Mammalian vector for expressing C-terminal 3xFlag tagged proteins.

Physical Map of pcDNA3.1-3xFlag:



| | |
|-------------|-----------------------------|
| Vector Name | pcDNA3.1-3xFlag |
| Vector Size | 5508 bp |
| Vector Type | Mammalian Expression Vector |

Rat Il2rb (XM_006241970) cDNA/ORF clone



Catalog Number: 708388-2

| | | | |
|-------------------|-------------------------------------|---------------------------------|------------|
| Expression Method | Constitutive, Stable / Transient | Antibiotic Resistance | Ampicillin |
| Promoter | CMV | Selection In Mammalian Cells | Neomycin |
| | | Protein Tag | Flag |