

Human DNAJB9 (NM_012328) cDNA/ORF clone



Catalog Number: 713396-3

General Information

Gene Name:

DnaJ heat shock protein family (Hsp40) member B9

Official Symbol: DNAJB9

Organism: Homo sapiens

RefSeq: NM_012328

Description

Sequence Description:

Identical with the Gene Bank Ref. ID sequence.

Vector: pCMV-SPORT6

Restriction Sites:

Shipping carrier:

Each tube contains approximately 5 µg - 10 µg of lyophilized plasmid.

Storage:

The lyophilized plasmid can be stored at ambient temperature for three months.

Quality control:

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

Sequencing primer list:

T7:TAATACGACTCACTATAGGG

M13 rev:CAGGAAACAGCTATGAC

Plasmid Resuspension protocol

1. Centrifuge at 5,000×g for 5 min.
2. Carefully open the tube and add 20 µl of sterile

water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than 5000×g.

5. Store the plasmid at -20 °C.

The plasmid is ready for:

Restriction enzyme digestion; PCR amplification; E. coli transformation; DNA sequencing


E.coli strains for transformation (recommended but not limited):

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5α and TOP10F'.

Vector Information

Vector for cloning and transient mammalian cell expression of cDNAs.

Physical Map of pCMV-SPORT6:



Vector Name

pCMV-SPORT6

Human DNAJB9 (NM_012328) cDNA/ORF clone



Catalog Number: 713396-3

Vector Size	4396 bp	Promoter	CMV
Vector Type	Mammalian Expression Vector	Antibiotic Resistance	Ampicillin
Expression Method	Transient	Selection In Mammalian Cells	-
		Protein Tag	None