

pLV3-CMV-KCNJ9(human)-3xFLAG-CopGFP-Puro

Name: pLV3-CMV-KCNJ9(human)-3xFLAG-CopGFP-Puro

Replication origin: ori

Vector ID: V024803

Fusion Tag: 3xFLAG


Antibiotic: Ampicillin

Expression method: Constitutive, Stable / Transient

Vector Type: Protein expression

Selection marker: Puro;CopGFP

Promoter: CMV



Plasmid Resuspension protocol

1. Centrifuge at 5,000×g for 5 min.
2. Carefully open the tube and add sterile water to dissolve the DNA: add 20 µl for 5 µg plasmid, and 100 µl for 100 µg plasmid.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than 5000×g.
5. Store the plasmid at -20 °C.