

pLV3-CMV-Esrrb(mouse)-6xHis-EGFP-EF1a-Puro

Name: pLV3-CMV-Esrrb(mouse)-6xHis-EGFP-EF1a-Puro

Vector ID: V032957

Antibiotic: Ampicillin

Vector Type: Protein expression

Promoter: CMV

Replication origin: ori

Primer 1: CMV-F


Fusion Tag: 6xHis-EGFP

Expression method: Constitutive, Stable / Transient

Selection marker: Puro;EGFP

Growth Temperature(°C): 37

Growth Strain(s): Stbl3



Plasmid Resuspension protocol

1. Centrifuge at 5,000×g for 5 min.
2. Carefully open the tube and add sterile water to dissolve the DNA: add 20 μl for 5 μg plasmid, and 100 μl for 100 μg plasmid.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than 5000×g.
5. Store the plasmid at -20 °C.