

HNRNPA1 antibody



Catalog Number: 111500

Product name

HNRNPA1 antibody

Immunogen

[Human HNRNPA1 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Pig; other species not tested.

Antibody description

HNRNPA1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HeLa cells, mouse colon tissue. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 34kd

Preparation

This antibody was obtained by immunization of HNRNPA1 recombinant protein (Accession Number: NM_002136). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111500(HNRNPA1 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:111500(HNRNPA1 antibody) at dilution of 1:100 (under 10x lens)