

Anti-VDAC1 Rabbit antibody



Catalog Number: 167500

Product name

Anti-VDAC1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

VDAC1 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human VDAC1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions


WB: 1/2000

IHC: 1/20-1/500


Validations



Western blot detection of VDAC1 in 293 cell lysates using VDAC1 Rabbit pAb(1:1000 diluted). Predicted band size:31KDa. Observed band size:31KDa.



Immunohistochemistry of VDAC1 in paraffin-embedded Human colon cancer tissue using



VDAC1 Rabbit pAb at dilution 1/50

Western blot detection of VDAC1 in Jurkat, Rat Brain, C6 cell lysates using VDAC1 Rabbit pAb(1:1000 diluted). Predicted band size:31kDa. Observed band size:31kDa.