

GRB2 antibody



Catalog Number: 111200

Product name

GRB2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GRB2 Rabbit Polyclonal antibody. Positive IF detected in HEK-293 cells. Positive WB detected in HEK-293 cells, A431 cells, C2C12 cells, rat brain tissue, rat spleen tissue. Observed molecular weight by Western-blot: 25-28kd

Preparation

This antibody was obtained by immunization of GRB2 recombinant protein (Accession Number: NM_002086). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100


Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111200(GRB2 antibody) at dilution of 1:1000



C2C12 cells (shcontrol and shRNA of GRB2) were subjected to SDS PAGE followed by western blot with Catalog No:111200(GRB2 antibody) at dilution of 1:1000. (Data provided by Angran Biotech (www.miRNAlab.com)).



Immunofluorescent analysis of HEK-293 cells

GRB2 antibody



Catalog Number: 111200

using Catalog No:111200(GRB2 Antibody) at

dilution of 1:25 and Alexa Fluor 488-conjugated
AffiniPure Goat Anti-Rabbit IgG(H+L)