

ARL3 antibody



Catalog Number: 108200

Product name

ARL3 antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

Specificity


Human, Mouse, Rat; other species not tested.

IF: 1:20-1:200

Antibody description

ARL3 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, human brain tissue, mouse brain tissue, mouse testis tissue, PC-3 cells, rat brain tissue. Positive IF detected in IMCD3 cells. Positive IHC detected in human testis tissue, human prostate cancer tissue. Observed molecular weight by Western-blot: 20 kDa

Validations



Preparation

This antibody was obtained by immunization of ARL3 recombinant protein (Accession Number: XM_017016260). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108200(ARL3 antibody) at dilution of 1:1000

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions

Recommended Dilution:



Immunohistochemical of paraffin-embedded human testis using Catalog No:108200(ARL3 antibody) at dilution of 1:50 (under 10x lens)




ARL3 antibody



Catalog Number: 108200

Immunohistochemical of paraffin-embedded human testis using Catalog No:108200(ARL3 antibody) at dilution of 1:50 (under 40x lens)



IF result of Catalog No:108200 (anti-ARL3,Green) in ciliated IMCD3 cells by Dr. Seongjin Seo. (Red is the primary cilium, and blue is the nucleus)