

Product name

CFP antibody

Specificity

Human; other species not tested.

Antibody description

CFP Rabbit Polyclonal antibody. Positive IHC detected in human spleen tissue. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 50 kDa

Preparation

This antibody was obtained by immunization of CFP recombinant protein (Accession Number: NM_002621). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

IHC: 1:20-1:200

Validations

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:109200(CFP Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human spleen tissue slide using Catalog No:109200(CFP Antibody) at dilution of 1:50 (under 40x lens)