

# Anti-AKT Rabbit antibody



Catalog Number: 168810

## Product name

Anti-AKT Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

AKT Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human AKT1/2/3

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, FC, IP

## Dilutions

WB: 1/2000


IHC: 1/50

ICC/IF: 1/20

FC: 1/50

IP: 1/20

## Validations



Western blot detection of AKT1/2/3 in Jurkat,Rat Brain,C6,CHO-K1,HeLa cell lysates using AKT Rabbit mAb(1:1000 diluted).Predicted band size:56kDa.Observed band size:56kDa.