

Anti-SCN9A antibody



Catalog Number: 175110

Product name

Anti-SCN9A antibody

Specificity

Human, Rat

Antibody description

Rabbit polyclonal antibody to SCN9A

Preparation

This antigen of this antibody was synthetic peptide within human nav1.7 aa 1950-1988 (cytoplasmic).

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 50% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ICC, IHC-P

Dilutions

ICC: 1:50-1:200

IHC-P: 1:50-1:200

Validations




Fig1: ICC staining NaV1.7 in SH-SY5Y cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Fig2: ICC staining NaV1.7 in A549 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Fig3: ICC staining NaV1.7 in HeLa cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.





Fig4: Immunohistochemical analysis of paraffin-embedded human placenta tissue using anti-Nav1.7 antibody. Counter stained with



hematoxylin.

Fig5: Immunohistochemical analysis of paraffin-embedded rat testis tissue using anti-Nav1.7 antibody. Counter stained with hematoxylin.