

# Human EPHA4/HEK8/SEK/TYRO1 recombinant protein



Catalog Number: 518510

## General Information

### Gene Name Synonym

Ephrin type-A receptor 4, HEK8, SEK, TYRO1, EPHA4, Tyrosine-protein kinase receptor SEK, Tyrosine-protein kinase TYRO1, EK8, hEK8, EPH-like kinase 8

### Protein Construction

A DNA sequence encoding the Human EPHA4/HEK8/SEK/TYRO1 Val20-Thr547 is expressed.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

### Stability

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

## Formulation

Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.

## Usage Guide


### Storage

Lyophilized protein should be stored at <-20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days.

### Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in ddH<sub>2</sub>O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

## SDS-PAGE



Human EPHA4/HEK8/SEK/TYRO1 recombinant protein