

DGKB antibody



Catalog Number: 109910

Product name

DGKB antibody

Immunogen

[Human DGKB Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

DGKB Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 91 kDa

Preparation

This antibody was obtained by immunization of DGKB recombinant protein (Accession Number: NM_145695). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109910(DGKB Antibody) at dilution of 1:1000