

# TRIM31 antibody



Catalog Number: 116310

## Product name

TRIM31 antibody

Recommended Dilution:

WB: 1:200-1:2000

## Immunogen

IP: 1:200-1:1000

[Human TRIM31 Recombinant protein \(His tag\)](#)

IHC: 1:20-1:200


## Specificity

Human; other species not tested.

## Validations

## Antibody description

TRIM31 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive IP detected in COLO 320 cells. Positive WB detected in COLO 320 cells. Observed molecular weight by Western-blot: 48 kDa



## Preparation

This antibody was obtained by immunization of TRIM31 recombinant protein (Accession Number: NM\_007028). Purification method: Antigen affinity purified.

COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:116310(TRIM31 antibody) at dilution of 1:400

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal


## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC, IP

## Dilutions



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:116310(TRIM31 antibody) at dilution of 1:50

# TRIM31 antibody



Catalog Number: 116310

IP Result of anti-TRIM31 (IP:Catalog No:116310, 4ug; Detection:Catalog No:116310 1:400) with COLO 320 cells lysate 3200ug.

