

GSTO1 antibody



Catalog Number: 107310

Product name

GSTO1 antibody

Immunogen

[Human GSTO1 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

GSTO1 Mouse Monoclonal antibody. Positive WB detected in Jurkat cells. Positive IHC detected in human liver tissue, human liver cancer tissue. Observed molecular weight by Western-blot: 27-30 kDa

Preparation

This antibody was obtained by immunization of GSTO1 recombinant protein (Accession Number: NM_004832). Purification method: Caprylic acid/ammonium sulfate precipitation.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107310(GSTO1 antibody) at dilution of 1:400



Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog No:107310(GSTO1 Antibody) at dilution of 1:50