

CIAPIN1 antibody



Catalog Number: 109310

Product name

CIAPIN1 antibody

Applications

ELISA, WB, IF

Immunogen

[Human CIAPIN1 Recombinant protein \(GST tag & His tag\)](#)

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

IF: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

CIAPIN1 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, HeLa cells, human brain tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 34kd

Preparation

This antibody was obtained by immunization of CIAPIN1 recombinant protein (Accession Number: NM_020313). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage


Store at -20°C. DO NOT ALIQUOT

Clonality


Polyclonal

Ig Type

Rabbit IgG



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:109310(CIAPIN1 antibody) at dilution of 1:400



Immunofluorescent analysis of HeLa cells using Catalog No:109310(CIAPIN1 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)