

# Anti-Cdk9 Rabbit antibody



Catalog Number: 167020

## Product name

Anti-Cdk9 Rabbit antibody

## Specificity

Human

## Antibody description

Cdk9 Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human Cdk9

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, IP

## Dilutions


WB: 1/1000

IHC: 1/20-1/100


ICC/IF: 1/20-1/50

IP: 1/20

## Validations



Immunocytochemistry of Cdk9(green) in HeLa cells using Cdk9 Rabbit pAb at dilution 1/50, and



Western blot detection of Cdk9 in K562,C6,3T3,HeLa cell lysates using Cdk9 Rabbit pAb(1:1000 diluted).Predicted band size:43kDa.Observed band size:43kDa.