

# ZNF671 antibody



Catalog Number: 117020

## Product name

ZNF671 antibody

## Immunogen

[Human ZNF671 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

ZNF671 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells. Observed molecular weight by Western-blot: 61-66 kDa

## Preparation

This antibody was obtained by immunization of ZNF671 recombinant protein (Accession Number: BC025728). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:117020(ZNF671 antibody) at dilution of 1:500