

TRUB1 antibody



Catalog Number: 116420

Product name

TRUB1 antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

Specificity

Human, Mouse; other species not tested.

Antibody description

TRUB1 Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue. Positive IP detected in mouse liver tissue. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 37 kDa

Preparation

This antibody was obtained by immunization of TRUB1 recombinant protein (Accession Number: NM_139169). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:1000


Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:116420(TRUB1 antibody) at dilution of 1:300



IP Result of anti-TRUB1 (IP:Catalog No:116420, 4ug; Detection:Catalog No:116420 1:300) with mouse liver tissue lysate 6400ug.




Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog

TRUB1 antibody

Catalog Number: 116420

No:116420(TRUB1 Antibody) at dilution of 1:50
(under 10x lens)



Immunohistochemistry of paraffin-embedded
human kidney tissue slide using Catalog
No:116420(TRUB1 Antibody) at dilution of 1:50
(under 40x lens)