

Anti-eEF1A1 Rabbit antibody



Catalog Number: 166820

Product name

Anti-eEF1A1 Rabbit antibody

WB: 1/5000

IHC: 1/20

Specificity

Human, Mouse, Rat

FC: 1/20

IP: 1/20

Antibody description

eEF1A1 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human eEF1A1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, FC, IP

Dilutions

Validations



Immunohistochemistry of eEF1A1 in paraffin-embedded Human tonsil using eEF1A1 Rabbit pAb



at dilution 1/20

Western blot detection of eEF1A1 in K562,C6,3T3,Hela cell lysates using eEF1A1 Rabbit pAb(1:1000 diluted).Predicted band size:50kDa.Observed band size:50kDa.