

ABLIM2 antibody



Catalog Number: 107820

Product name

ABLIM2 antibody

Immunogen

[Human ABLIM2 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ABLIM2 Rabbit Polyclonal antibody. Positive IHC detected in human skeletal muscle tissue, human kidney tissue. Positive WB detected in mouse brain tissue, mouse skeletal muscle tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 60-70 kDa

Preparation

This antibody was obtained by immunization of ABLIM2 recombinant protein (Accession Number: NM_032432). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


IP: 1:500-1:5000

IHC: 1:20-1:200

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107820 (ABLIM2 antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human skeletal muscle slide using Catalog No:107820 (ABLIM2 Antibody) at dilution of 1:50

ABLIM2 antibody



Catalog Number: 107820

IP Result of anti-ABLIM2 (IP:Catalog No:107820, 4ug; Detection:Catalog No:107820 1:1000) with mouse brain tissue lysate 5200ug.

