

Anti-SMN/Gemin 1 Rabbit antibody



Catalog Number: 168820

Product name

WB, ICC/IF

Anti-SMN/Gemin 1 Rabbit antibody

Dilutions

Specificity

1:1000-1:5000

Human

1:20-1:50

Antibody description

SMN/Gemin 1 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Gemin 1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality


Monoclonal

Ig Type


Rabbit IgG

Applications

Validations



Western blot detection of SMN/Gemin 1 in K562, HeLa cell lysates using SMN/Gemin 1 Rabbit mAb(1:1000 diluted). Predicted band size: 32kDa. Observed band size: 35kDa.



Immunocytochemistry of Gemin 1 (green) in HeLa cells using Gemin 1 Rabbit mAb at dilution 1/200, and DAPI (blue)