

Phospho-TBC1D4-T642 pAb



Catalog Number: 166520

Product name

Phospho-TBC1D4-T642 pAb

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to Phospho-TBC1D4-T642

Preparation

Antigen: A synthetic phosphorylated peptide around T642 of human TBC1D4 (NP_055647.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - Phospho-TBC1D4-T642 pAb

Western blot analysis of extracts of 293T cells, using Phospho-TBC1D4-T642 antibody at 1:2000 dilution. 293T cells were treated by IGF-1 (50ng/ml) for 5 minutes after serum-starvation overnight. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% BSA. Detection: ECL Enhanced Kit. Exposure time: 30s.