

**Product name**

LIAS antibody

**Immunogen**

[Human LIAS Recombinant protein \(GST tag & His tag\)](#)

**Specificity**

Human, Mouse, Rat; other species not tested.

**Antibody description**

LIAS Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 42kd

**Preparation**

This antibody was obtained by immunization of LIAS recombinant protein (Accession Number: NM\_006859). Purification method: Antigen affinity purified.

**Formulation**

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

**Storage**

Store at -20°C. DO NOT ALIQUOT

**Clonality**

Polyclonal

**Ig Type**

Rabbit IgG

**Applications**


ELISA, WB, IHC

**Dilutions**


Recommended Dilution:

WB: 1:500-1:5000

IHC: 1:20-1:200

**Validations**

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112220(LIAS antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:112220(LIAS antibody) at dilution of 1:100 (under 40x lens)